



TECHNICAL SPECIFICATION



ACCENT BUDGET RANGE BED

Product Features:

- ✓ Fully Electrically Operated
- ✓ Reduced Shear and Friction
- ✓ Reliable and Safe
- ✓ Wooden Integral Siderails
- ✓ Light Beech Finish
- ✓ Metal Mattress Platform
- ✓ Simple Handset Controls
- ✓ Meets Industry Safety Standards

VALUE

SIDERAILS
INCL.

Technical:

Safe Working Load	180kg
Outside Width	1020mm
Inside Width	900mm
Height Range	400 - 800mm
Total Weight	64kg
Backrest Angle	0 - 70 degrees
Length Inside	2240mm
Length Outside	2080mm

Warranty:

Frame	N/A
-------	-----

FOR MORE INFORMATION OR TO ORDER CALL 01423 331000 NOW