

Burn Gel Dressings and Tubes

Description

AeroBurn™ is a market leading burn relief range for the immediate first aid treatment of burns, scalds and sunburn. AeroBurn™ Gel offers immediate and rapid cooling to help prevent further tissue damage, relieving pain and aiding in the healing process.

AeroBurn™ 3.5g satchets offer a sterile single use application while 25 and 50g tubes are ideal for large application, multiple uses and redressing. AeroBurn™ Gel Infused dressings absorb and dissipate the heat associated with burns, reduce inflammation and stop the burning to set you on your way to healing.

Features & Benefits

Enriched with Aloe Vera

- The AeroBurn™ Gel formula is enriched with the natural healing powers of Aloe Vera. Aloe Vera speeds up the recovery of skin wounds and is proven to aid cellular rejuvenation.

Contains Melaleuca Oil

- Melaleuca Oil, a natural product of the tea tree, is a proven antiseptic and antifungal remedy which is beneficial for all skin types and does not alter skin cell rejuvenation. Melaleuca oil also has local anaesthetic and anti-inflammatory properties.

Non-Woven Carrier

- AeroBurn Gel infused Dressing use a Non Woven Carrier, which has superior thermodynamics when compared to any other foam carrier. Results show that the AeroBurn Non-Woven carrier dressing absorbs 50% more than an equivalent sized foam carrier, helping retain more gel at the wound site to cool the burn and soothe pain.

Non-Adherent, Non-Toxic


- AeroBurn Burn Gel and Dressings are non-adherent and non-toxic, ensuring minimal irritation and maximum relief for the wound.

Easy to Use

- AeroBurn Burn Gel Dressings offer more effective relief for burn wounds, as well as being more cost-effective. The burn dressing can simply remain in situ for up to 2 hours, providing relief with minimal attention, redressing or cooling using running water.

Contraindications/Precautions

- This product is not a substitute for professional medical attention from a healthcare practitioner
- Medical attention should be sought for any burns, larger than a 20 cent piece; involving airways, eyes, hands, face or genitals; associated with inhaled smoke, gas, fumes or caused by electricity; deeper than superficial or where severity is uncertain.



Code	Item	Size Description	Items/Unit	MOQ	Shipper
AB3	AeroBurn Burn Gel Satchet	3.5g	satchet/1	500	1000
ABD10	AeroBurn Burn Dressing	10cm x 10cm	dressing/1	36	144
ABS50	AeroBurn Burn Spray	50ml	spray/1	6	96