



TECHNICAL SPECIFICATION



AGENA BARIATRIC PROFILING BED

Product Features:

- ✓ Electrically operated backrest, height, tilt and knee break
- ✓ 3 Detachable Aluminium side rails
- ✓ 4 lockable twin castors
- ✓ Easy to clean and low maintenance
- ✓ Weight Limit 318kg (50 stone)
- ✓ 27 - 68cm height variation
- ✓ Trendelenburg function as standard
- ✓ Remote handset control
- ✓ 27cm from floor to top of platform



Technical:

| | |
|--------------------------|-----------------------|
| <u>Safe Working Load</u> | <u>318kg</u> |
| <u>Outside Width</u> | <u>1360mm</u> |
| <u>Inside Width</u> | <u>1200mm</u> |
| <u>Height Range</u> | <u>270 -680mm</u> |
| <u>Total Weight</u> | <u>N/A</u> |
| <u>Backrest Angle</u> | <u>0 - 72 degrees</u> |
| <u>Length Inside</u> | <u>1980mm</u> |
| <u>Length Outside</u> | <u>2180mm</u> |

Warranty:

| | |
|-----------------|---------------|
| <u>Actuator</u> | <u>1 Year</u> |
| <u>Frame</u> | <u>1 Year</u> |

FOR MORE INFORMATION OR TO ORDER CALL 01423 331000 NOW