



TECHNICAL SPECIFICATION



CENTAURI EXTRA LOW PROFILING BED

Product Features:

- ✓ Drops to Just 180mm from the Floor
- ✓ Four Section Profiling
- ✓ Trendelenburg Function
- ✓ Four Braked Castors
- ✓ Does Not Require Siderails
- ✓ CE Marked
- ✓ Metal Mesh Mattress Platform
- ✓ Available in Cherry
- ✓ Enclosed Actuators
- ✓ Simple Handset Controls

LOW
BED

NO HAND
RAILS

Technical:

Safe Working Load	180kg
Outside Width	1040mm
Inside Width	900mm
Height Range	180 - 560mm
Total Weight	70kg
Backrest Angle	0 - 75 degrees
Length Inside	2030mm
Length Outside	2130mm

Warranty:

Actuator	3 Years
Frame	3 Years

FOR MORE INFORMATION OR TO ORDER CALL 01423 331000 NOW