



andway
healthcare in hand

Glove Size Guide

To maintain protection from contamination, dexterity and sensitivity for your hands.

Please note you should use the diagram below as a guide but should try on the gloves to confirm you have a good fit.



PALM WIDTH	EXAMINATION GLOVE SIZE
up to 80 mm	Extra Small (XS)
70 - 90 mm	Small (S)
85 - 105 mm	Medium (M)
100 - 120 mm	Large (L)
≥ 120 mm	Extra Large (XL)

PALM WIDTH	PROCEDURE GLOVE SIZE
63 - 71 mm	5
68 - 76 mm	5.5
72 - 83 mm	6
78 - 88 mm	6.5
84 - 94 mm	7
90 - 100 mm	7.5
94 - 108 mm	8
102 - 114 mm	8.5
108 - 120 mm	9
115 - 127 mm	9.5

To view our full range of Medical Gloves, visit andwayhealthcare.com

T: 0333 600 1400

E: sales@andwayhealthcare.org.uk