



## Product Specifications Data Sheet

<b>Product Name</b>	Examination Gloves (Medical)		
<b>Product Description</b>	Blue Soft Vinyl PF Gloves		
<b>Range</b>	Premier		
<b>Manufacturer</b>	Andway Healthcare		
<b>Marking</b>	Thin lightweight design has a natural feel for wearers and patients/clients alike. Easy to don and comfortable to wear for longer periods. Smooth surface creates minimal friction therefore ideal for close personal contact in health care and beauty treatment where there is risk of blood or bodily fluid contact.		
<b>Material</b>	Vinyl		
<b>Sterilisation</b>	Non-sterile		
<b>Issue Date</b>	20 July 2023	<b>Version No</b>	1



<b>Primary use</b>	A single use device, intended to be used as a protective barrier when worn on the hands of healthcare providers during patient examinations/treatment.
--------------------	--

Packaging									
Product code	Size	Palm Width (mm)	Colour	Length (mm)	Palm Thickness (mm)	Force at Break	Items Per Pack	Qty- Inner	Qty – Shipper Carton
VPFBG/SM	Small	85±5	Clear	240	0.07	3.6N (min)	10	100	1000
VPFBG/MED	Medium	95 ± 5							
VPFBG/LGE	Large	105 ± 5							
VPFBG/XL	Extra Large	115 ± 5							

<b>Quality certification</b>	ISO 13485
<b>Device Classification</b>	Class I – Medical Device Directive 93/42/EEC
<b>Conformity to Guidelines</b>	EN 455 Medical Gloves for single use. Parts 1-4 AQL 1,5
<b>Storage conditions</b>	Store in a cool dry place. Keep away from heat and direct sunlight
<b>Shelf life</b>	5 years from the date of manufacture with the above storage condition
<b>Additional information</b>	Smooth Surface. Ambidextrous Fit.

*Medical Device and Reporting System:*

*According to the official criteria of MDD, incidents caused by this device must be reported immediately to our offices.*

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

Andway Healthcare Ltd	Atlantis House, Manse Lane
T: 0333 600 1400	Knarborough
W:www.andwayhealthcare.com	HG5 8LF