

PRODUCT SPECIFICATION

Product name	Latex Examination Gloves
Range	Premier Performer
Type	Powdered, Non-Sterile
Marking	Product name, Premier Logo, Size, Code Ref, Lot Number, Expiry date, Latex Free symbol, non-sterile symbol, CE Mark, and single use symbols
Shape	Standard cuff length and ambidextrous.
Material	Latex
Colour	Natural



Primary use	A single use, non-sterile examination glove intended to be used for medical purposes to prevent contamination between patient and examiner.
--------------------	---

Ref	Size	Length (mm)	Palm width (mm)	Thickness (mm) (minimum at palm)	Force at Break (N)
P7401	S	≥240	80 ± 10	0.10-0.11	≥ 6.0 (during shelf life) ≥ 6.0 (after challenge testing)
P7402	M	≥240	95 ± 10	0.10-0.11	≥ 6.0 (during shelf life) ≥ 6.0 (after challenge testing)
P7403	L	≥240	110 ± 10	0.10-0.11	≥ 6.0 (during shelf life) ≥ 6.0 (after challenge testing)
P7404	XL	≥240	≥ 110	0.10-0.11	≥ 6.0 (during shelf life) ≥ 6.0 (after challenge testing)

Specifications according to EN 455

Packaging	Dispenser box S; M; L; XL.	100 Gloves per dispenser
	Outer Carton	1000 (10 Dispenser Boxes of 100 Gloves)

Quality certification	ISO 13485
Conformity to Regulations	Class I – European Medical Device Regulation 2017/745
Conformity to standards	EN 455 Medical Gloves for single use Parts 1-4 AQL 1,5
Storage conditions	Store in a cool dry place away from heat and direct sunlight
Shelf life	Refer to product packaging for date of expiry
Additional information	Finger textured to provide excellent reliability

Medical Device and Reporting System:

Incidents caused by this device must be reported immediately to our offices.

Andway Healthcare
Atlanits House, Manse Lane,
Knaresborough, HG5 8LF

T: 0333 600 1400
E: sales@andwayhealthcare.org.uk
W: andwayhealthcare.com